


PC-based Speech Audiometer



Oscilla[®] Speech



Full featured
Easy to use
Ultra portable
Single USB connection



The new Oscilla USB-350SP speech audiometer combines a full-featured diagnostic PC-based audiometer with speech, ease of use and ultra portability with a single USB connection.

With 24-bit digital audio quality and the AudioConsole software, the USB-350SP masters live-voice and WAV-file test with several masking types such as white noise, speech noise and narrow band noise.

The USB-350SP is fully customizable to suit your individual needs and is built on a future proof technical platform making it ideal for adaption / extension.

Includes:

Patient headset, operator headset, integrated patient response, bone conductor, AudioConsole software, user manuals, carrying bag



Manufactured by Inmedico A/S Denmark
Møgelhøj 2 - DK-8520 Lystrup

www.inmedico.com - info@inmedico.com
tel.: +45 86 74 26 22 - fax: +45 86 74 26 16

